

DABIRAIR™ CONTROLLERS


DABIR SURFACES: ADVANCED HUMAN INTERFACE TECHNOLOGIES

dabir-surfaces.com

CONTROLLERS



MODEL CA-1001



DabirAIR Controllers are engineered to address both caregiver and patient needs by incorporating advanced electronics and proprietary sensing technology to control, monitor and communicate. Like all Dabir products, these units are intuitive, simple to operate and extremely reliable. Key system features include wireless communication between overlay and controller, patient usage-profile monitoring and reporting and solid-state touch input keys with illuminated graphics, all in a quiet, lightweight and easy to store design.


- Easy Power Activation
- Smart Interconnect System
- Start/Stop Touch Operation
- Variable Cycle Speed Settings
- Variable Firmness Settings

ONE MISSION:
TO REDUCE THE RISK OF PRESSURE ULCERS

 **DABIR™**
SURFACES

SPECIFICATIONS

MODEL	CA-1001
Supply Voltage	100-240 VAC
Supply Frequency	50-60 Hz
Power Input	170VA Max
Size	356 x 178 x 127mm (14 x 7 x 5 in.)
Weight	3.95 kg (8.7 lb)
Type of Protection Against Electric Shock	Class 1
Degree of Protection from Electrical Shock	Type BF
Environment Conditions:	
– Operation: 10°C to 35°C (50°F to 95°F) 30-80% RH	
– Storage / Transport: -40°C to 60°C (-40°F to 140°F) 10-95% RH	
– Maximum Operating Altitude: 3000m (9842 ft.)	



No latex is used in the manufacturing of this product
Data subject to change



Follow instructions for use

DABIRAIR™ ACCESSORIES

MODEL	DESCRIPTION
CA-9001	Replacement Hose Assembly – 9'
CA-9002	Replacement Power Cord – 10'
Pending	IV Pole Mounting Kit
Pending	MedSurg Bed Rail Mounting Kit

DabirAIR Accessories are designed to address the ever changing needs of patient care, including various options for controller mounting, optimal cord / hose routing solutions and hardware replacement.

We encourage customer feedback so that we may continuously improve our products and services. Working together, we can achieve our goal of reducing pressure ulcer risk and improving patient outcomes.

24585 Evergreen Road, Southfield, MI 48075 USA
Tel: +1 (248) 796-0802
Fax: +1 (248) 354-6667

SALES

sales@dabir-surfaces.com

TECHNICAL AND WARRANTY SUPPORT

support@dabir-surfaces.com
Tel: +1 (885) 478-0694



dabir-surfaces.com