


DABIRAIR™

MICROCLIMATE

DABIR SURFACES: ADVANCED HUMAN INTERFACE TECHNOLOGIES

dabir-surfaces.com

MICROCLIMATE



MODEL MC-303584 MEDSURG OVERLAY



DabirAIR Microclimate Overlays combine the reduced skin shear and tissue relieving properties of Dabir's patent pending alternating Micropressure Overlay with a low-air-loss microclimate feature. This enhancement further reduces the risk of tissue breakdown due to excess heat and moisture which can lead to skin maceration and increased friction during immobilized patient turning.

- Multi-Patient / Semi-Disposable
- Alternating Nodal Support
- Reduced Skin Shear
- Removes Excess Heat/Moisture
- Easy Wipe-Down Cleaning
- Covers with Standard Sheets

ONE MISSION:
TO REDUCE THE RISK OF PRESSURE ULCERS

 **DABIR™**
SURFACES


SPECIFICATIONS

MODEL	MC-303584
Cleanable	Yes (Standard Wipes)
Disposable (Optional)	Yes
Self-Contouring	Yes
Microclimate Feature	Optional
Usage Life	1,440 Hrs. (Cumulative)
Typical Replacement	60 Days (In Application)
Length	2134 mm (84 in.)
Width	838 mm (33 in.)
Height (Deflated)	2 mm (0.1 in.)
Height (Inflated)	25 mm (1 in.)
Weight	2.95 kg (6.5 lbs.)
Intelligent / Wireless Data Collection	RFID Technology
Warranty	1 Year Limited


No latex is used in the manufacturing of this product
Data subject to change




Follow instructions for use



Simple ergonomic connection



Covers with standard bed sheet
(Customer supplied. Not shown)



We encourage customer feedback so that we may continuously improve our products and services. Working together, we can achieve our goal of reducing pressure ulcer risk and improving patient outcomes.



24585 Evergreen Road, Southfield, MI 48075 USA
Tel: +1 (248) 796-0802 Fax: +1 (248) 354-6667
dabir-surfaces.com

SALES AND CUSTOMER SUPPORT

Tel: +1 (888) 559-3641
sales@dabir-surfaces.com

TECHNICAL AND WARRANTY SUPPORT

Tel: +1 (888) 559-3642
support@dabir-surfaces.com